



binding of SKOV-3 particles to the surface-bound LXY30 ligands (Figure 2, Table 1). In this study, only three concentrations of analyte were measured, however, it is generally recommended to have five or more concentrations to evaluate affinity and kinetics.

Thickness and refractive index (RI) of the adsorbed EV layer were evaluated by fitting the full SPR curves in the LayerSolver™ software (Figure 3). Analysis was based on multi-wavelength measurement of MP-SPR, provided by 670 and 758 nm laser sources. The surface mass density of surface bound EVs was estimated to be 90 ng/cm<sup>2</sup>. Refractive index of the EV layer was 1.341, which is slightly higher than the RI of aqueous buffer solution like PBS. The calculated thickness of the exosome layer was only 34 nm based on MP-SPR, even if average diameter of the particles was 177 ± 87 nm based on NTA. This might be explained by the fact that MP-SPR provides average thickness from the measurement area, while the EVs were unevenly distributed over the sensor. The low thickness value was most likely caused by empty spaces between the exosomes bound onto the surface. Also, the deformation of the EVs shape due to binding may explain the observed lower value of thickness. Furthermore, size and concentration of the extracellular vesicles can be assessed with high accuracy when using the obtained MP-SPR data with a modelling protocol published by Rupert *et al.* 2016.

## Conclusions

MP-SPR has proven to be an excellent platform to characterize interactions of extracellular vesicles, as shown in this example. A single MP-SPR measurement provided several numerical values, such as affinity and kinetic constants, thickness, refractive index and adsorbed mass, while the experiment was performed in a label-free manner. The approach presented here is not limited to extracellular vesicles but is applicable to other types of nanoparticles as well, including polymeric and metallic nanoparticles. The possibility to work with various samples, surfaces and matrices make MP-SPR an outstanding platform for biosensor development and life science research.

See also how MP-SPR can measure the formation of protein corona on liposome in 100% serum Application Note #151 and how exosomes and other nanoparticles uptake by living cells was measured Application Note #156.

### Original publication:

Carney *et al.* Advanced Biosystems 2017, Volume 1, Issue 5, 1600038

### Reference:


Rupert *et al.* Analytical Chemistry, 2016, Vol. 88 (20), p. 9980–9988

### Recommended instrumentation for reference assay experiments

MP-SPR Navi™ 200 OTSO, 400 KONTIO, 210A VASA, 410A KAURIS, 220A NAALI and 420A ILVES with additional wavelengths

Sensor surfaces: CMD, Au or inorganic coating


Software: MP-SPR Navi™ Control, DataViewer, LayerSolver and TraceDrawer for MP-SPR Navi™



**Figure 2.** Affinity and kinetic constants of peptide – extracellular vesicle interaction was determined using TraceDrawer™ for MP-SPR Navi™. The solid black line shows measured data, and the dashed orange line represents the fit to sensogram from sequential injections of three different populations of EVs having A) 1%, B) 5%, and C) 10% coverage of α3β1 integrin.

Estimated α3β1	$k_{a1}$ (1/M <sup>2</sup> s)	$k_{d1}$ (1/s)	$K_{D1}$ (M)	$k_{a2}$ (1/M <sup>2</sup> s)	$k_{d2}$ (1/s)	$K_{D2}$ (M)
1%	$5.14 \times 10^5$ (±6.97×10 <sup>3</sup> )	$9.35 \times 10^{-3}$ (±9.11×10 <sup>-5</sup> )	$1.82 \times 10^{-8}$ (±4.23×10 <sup>-10</sup> )	$1.28 \times 10^6$ (±7.28×10 <sup>2</sup> )	$7.00 \times 10^{-5}$ (±1.88×10 <sup>-5</sup> )	$5.46 \times 10^{-11}$ (±1.47×10 <sup>-11</sup> )
5%	$6.32 \times 10^5$ (±1.85×10 <sup>4</sup> )	$7.30 \times 10^{-3}$ (±1.54×10 <sup>-4</sup> )	$1.16 \times 10^{-8}$ (±5.83×10 <sup>-10</sup> )	$2.52 \times 10^5$ (±5.93×10 <sup>2</sup> )	$6.05 \times 10^{-5}$ (±2.53×10 <sup>-5</sup> )	$2.40 \times 10^{-10}$ (±1.01×10 <sup>-10</sup> )
10%	$2.90 \times 10^5$ (±1.86×10 <sup>4</sup> )	$7.40 \times 10^{-3}$ (±1.26×10 <sup>-4</sup> )	$2.55 \times 10^{-8}$ (±2.08×10 <sup>-9</sup> )	$1.26 \times 10^5$ (±6.30×10 <sup>2</sup> )	$6.07 \times 10^{-5}$ (±2.13×10 <sup>-5</sup> )	$4.82 \times 10^{-10}$ (±1.71×10 <sup>-10</sup> )

**Table 1.** Affinity and kinetic constants of SKOV-3 particles binding to LXY30 peptide.



**Figure 3.** Full MP-SPR curves before and after adsorption of SKOV-3 extracellular vesicles were analyzed using LayerSolver software to solve thickness and refractive index of the particle layer. Multi-wavelength analysis was performed but for presentation clarity only 670 nm measurement curves are presented here.